

DAQ System with BioWare® 4.0

Type 5691A1

Data Acquisition and Analysis System for Biomechanics

Data acquisition system for connecting and controlling two multicomponent force plates with integral charge amplifiers. The system is connected to a USB 2.0 port of the PC and operated with the included software BioWare 4.0.

- Easy of installation with USB 2.0
- Remote control of integral charge amplifiers
- Powerful data acquisition and signal processing
- Versatile data analysis and filters
- External and internal trigger

Description

The DAQ system with BioWare 4.0 consists of a connecting box for one or two Kistler multicomponent force plates and one integral 16-bit A/D converter to digitize the plates' analog output signals. The system is connected to a USB 2.0 port of the PC. The integral charge amplifiers of the connected Kistler force plates are supplied via the connecting box and controlled by means of the supplied software (measuring range and reset/operate).

Application

The Type 5691A1 with BioWare 4.0 is designed specifically to fully exploit the capabilities of Kistler's piezoelectric force plates in biomechanics applications. The 16-bit resolution of the measurement signals and high sampling rate of up to 17 kS/s in conjunction with Kistler force plates allow a very wide range of applications. The system as a whole is therefore equally ideal for measuring highly dynamic processes, very small measurands and slow phenomena. The additional options of acquiring any analog signals rather than just those from force plates, and external triggering with pre- and post-trigger capability, underscore the versatility of the system for use in basic research, sports science, gait analysis, neurology, ergonomics, etc., etc.

Technical Data

General Data

| | | |
|-----------------------------|----|------------|
| Dimensions | mm | 208x65x250 |
| Total weight | kg | 2,05 |
| Operating temperature range | °C | 0 ... 50 |

Power Supply Voltage

| | | |
|-------------------|-----|-----------|
| Power supply | VDC | 11 ... 15 |
| Power consumption | VA | 6 |



A/D-Converter

| | | |
|--|--------------------|-----------------|
| Number of channels | | 16 |
| Resolution (per channel) | Bit | 16 |
| Input voltage range (software selectable) | V | ±1, ±2, ±5, ±10 |
| Sampling rate (software selectable) | S/s | 0,6 ... 50 000 |
| | max. @ 2 channels | kS/s 50 |
| 1 Force plate | max. @ 8 channels | kS/s 17 |
| 2 Force plates | max. @ 16 channels | kS/s 9,5 |

Connections

USB 2.0

| | |
|----------------------------|--------------------|
| USB In (uplink, to the PC) | USB Type B, female |
| USB Out (downlink, free) | USB Type A, female |

Force Plate 1/2

| | | |
|----------------------|---|-----|
| Input voltage (max.) | V | ±15 |
|----------------------|---|-----|

External Trigger

| | | |
|----------------------------------|---------------------|--------------|
| Input voltage | | BNC neg. |
| Pull-Up resistance 10 kΩ to ±5 V | | |
| | max. | VDC 12 |
| | high or input open | VDC >3,6 |
| | low | VDC <0,6 |
| Trigger mode | standard | rising edge |
| | software selectable | falling edge |

Conforms to the **CE** safety standards (73/23/EG) for electrical equipment and systems:

EN 60601-1:2005, EN 61010-1:2001
and the EMC standards (89/336/EG):
EN 60601-1:2005 (EN 55022 Class B), EN 61000-6-3:2004
(EN 55022 Class B), EN 61000-6-4:2001 (EN 55011 Class B),
EN 60601-1:2005, EN 61000-6-1:2001, EN 61000-6-2:2005

BioWare®

BioWare software is the engine behind the force plate system. It collects data from the force plates, converts the trials into useful information and plots the results. The force plates and charge amplifiers are fully remote controlled by BioWare thus making the system extremely flexible and easy-to-use.

Parameters of Gait




Fig. 1: Ground reaction forces (GRF)




Fig. 2: Center of pressure (COP)




Fig. 3: Frictional torque T_z




Fig. 4: Force vector

Other functions

- Coefficient of friction (COF)
- Frequency analysis, statistics, digital filters
- Full Windows® functionality

BioWare provides several performance specific evaluations.

Parameters of Countermovement Jump CMJ




Fig. 5: Jump force




Fig. 6: Power




Fig. 7: Jump height (COM)




Fig. 8: Force gradient (Explosivity)

Other parameters

- Acceleration, velocity and displacement of the center of mass (COM)
- Work, energy, impulse
- Statistics, digital filters

5691A_000-633e-09.08

Dimensions





Fig. 9: DAQ system with BioWare Type 5691A

Pin Allocation D-Sub37, male (Force Plate 1/2)

| | | | | |
|----|--------------------------|----|--------------------------|----------------------|
| 1 | Exct. +12 VDC | 20 | Data IO5 (reserve) | |
| 2 | n.c. | 21 | Data IO6 (reserve) | |
| 3 | n.c. | 22 | B Range Select Group I | |
| 4 | A Range Select Group I | 23 | Operate/NotReset | |
| 5 | n.c. | 24 | Data IO7 (reserve) | |
| 6 | n.c. | 25 | n.c. | |
| 7 | Exct. GND | 26 | n.c. | |
| 8 | n.c. | 27 | B' Range Select Group II | |
| 9 | A' Range Select Group II | 28 | Control GND | |
| 10 | n.c. | 29 | Control GND | |
| 11 | Signal GND | - | Force Plate 1 | Force Plate 2 |
| 12 | Signal GND | 30 | CH8 (Fz4) | CH16 (Fz4) |
| 13 | Signal GND | 31 | CH7 (Fz3) | CH15 (Fz3) |
| 14 | Signal GND | 32 | CH6 (Fz2) | CH14 (Fz2) |
| 15 | Signal GND | 33 | CH5 (Fz1) | CH13 (Fz1) |
| 16 | Signal GND | 34 | CH4 (Fy23) | CH12 (Fy23) |
| 17 | Signal GND | 35 | CH3 (Fy14) | CH11 (Fy14) |
| 18 | Signal GND | 36 | CH2 (Fx34) | CH10 (Fx34) |
| 19 | Signal GND | 37 | CH1 (Fx12) | CH9 (Fx12) |



Typical Measuring Chain




Fig. 10: Configuration of a typical measuring chain with Kistler DAQ system with BioWare®

5691A_000-633e-09_08

DAQ-System with BioWare Version 4.0 Type 5691A1

Data Acquisition and Analysis Tool for Biomechanics

USB 2.0, for max. 2 Force Plates, with built-in charge amplifier

- USB 2.0 DAQ system with BioWare Type 5691A (16 channels, 16 bit)
- BioWare software

System Requirements

- Microsoft Windows® 2000 (SP3), or Windows XP operating system
- Intel Pentium III class processor (500 MHz or higher recommended)
- 512 MB of RAM
- Video Display set to at least 800x600, 256 colors, small fonts selectet
- Min 125 MB of free disk space
- Microsoft compatible mouse
- Windows Installer Version 1.1 or later
- Adobe® Acrobat® Reader®
- 1 free USB 2.0 port

Included Accessories

- | | Type/Art. No. |
|--|---------------|
| • USB 2.0 connecting cable, length 2 m | – |
| • Universal AC/DC adapter, 100 ... 240 V~ 12 VDC | 5.510.276 |
| • Self-adhesive base, black, 20,5x7,6 mm | 5.211.368 |
| • BioWare software CD-ROM | 2812A-04-2 |
| • HASP licence key (USB) | – |
| • Instruction manual | 2812A_002-312 |

Optional Accessories

- | | Type/Art. No. |
|---|---------------|
| • Connection cable for | |
| – Force platforms w/ integr. amplifier (straight connector) | 1758A... |
| – Force platforms w/ integr. amplifier (angle connector) | 1759A... |
| – External charge amplifier Type 9865E... | 1769A1 |
| – External control unit Type 5233A2 | 1500B5 |
| – Analog signals (8x BNC pos.) | 1500A67 |
| • Mounting kit consisting of 2 brackets and 4 screws | 7.511.339 |

Ordering Code

- | | |
|-------------------------------|-------------|
| • DAQ system with BioWare 4.0 | Type 5691A1 |
|-------------------------------|-------------|